

A TSKgel ODS-100Z 3 μ m (4.6 mm I.D. \times 25 cm) analytical column was used. ADAM (9-anthryldiazomethane) was used as the fluorescent derivatization reagent, and after reacting at 40°C, separation was performed.

Figure 2 shows a chromatogram of a standard sample. Calibration curves having good linearity in the concentration range 0.5 to 50 mg/L were obtained for both EPA and DHA. Under these analytical conditions, the limit of quantitation was 0.02 mg/L. Figure 5 shows a chromatogram from an analysis of EPA and DHA within a fish-oil-containing supplement. The quantity of each analyte was determined without any effects of contaminants, and the recovery was 92 to 95%.

- 1) 10 mg of ADAM was dissolved in 10 mL of ethyl acetate
↓
- 2) 10 mL of ethyl acetate was added to 20 mg of sample and dissolved with sonicator
The sample solution was filtered (PTFE, 0.5 μ m pore)
The filtrate was diluted 10-fold with ethyl acetate
↓
- 3) 200 μ L of 1) and equal volume of 2) was mixed and reacted at 40°C for 60 min
↓
- 4) reaction solution was diluted 10-fold with acetonitrile and applied to HPLC

Figure 3 Procedure of pretreatment





Figure 4 Calibration curves




Figure 5 Chromatogram of fish-oil-containing supplement