

# DEC Caps for use with Gilson ASPEC™ Automated SPE Systems

- 1 and 3ml options
- Compatible with all Gilson ASPEC™ Systems
- Excellent sealing – needle and SPE Column
- Compatible with all 1 and 3ml Syringe Barrel SPE Columns
- Suitable for use with tab-less and conventional or winged/tabbed SPE columns



## Introduction

Kinesis DEC Caps allow the automation of SPE methods using Gilson ASPEC™ Systems. Fully compatible with the majority of SPE columns, the Kinesis DEC Caps are inserted into the inlet of the syringe barrel style SPE column, allowing a positive pressure SPE method to be utilised.

## 1 and 3ml Options

Kinesis offers both 1 and 3ml DEC Caps. With the majority of analytical SPE methods using 1 or 3ml SPE Columns, the Kinesis DEC Cap supports most clinical, environmental, agrochemical and food solid phase extraction methods.

## Tab-less & Conventional/Tabbed/Winged Columns

The outer diameter allows the Kinesis DEC Caps to be used with either conventional/tabbed/winged SPE Columns, as well as the tab-less version, often used in 96- and 48-well plate footprints.

## Optimised Material

The proprietary material and design ensures an excellent seal with both the internal surface of the SPE Column and the needle of the appropriate Gilson ASPEC System.

## Ordering Information

| Part Number | Description                            | Pack Size |
|-------------|----------------------------------------|-----------|
| DC-GIL-001  | Kinesis DEC Caps for Gilson ASPEC, 1ml | 1000      |
| DC-GIL-003  | Kinesis DEC Caps for Gilson ASPEC, 3ml | 1000      |

## Samples

Samples of Kinesis 1 and 3ml DEC Caps are readily available. Please contact your local Kinesis office or distributor.



### UNITED KINGDOM

#### Kinesis Ltd

Tel: +44 (0)1480 212122  
Fax: +44 (0)1480 212111  
sales@kinesis.co.uk  
Web: kinesis.co.uk

### GERMANY, SWITZERLAND & AUSTRIA

#### Kinesis GmbH

Tel: +49 (0)2173 89 05-0  
Fax: +49 (0)2173 89 05-77  
sales@kinesisgmbh.de  
Web: kinesisgmbh.de

### USA & CANADA

#### Kinesis Inc

Tel: (518) 289-5817  
Toll free in USA: (866) 934-6353  
Fax: (518) 289-5818  
sales@kinesis-usa.com  
Web: kinesis-usa.com

### AUSTRALIA & NEW ZEALAND

#### Kinesis Australia Pty Ltd

Tel: +61 (0)7 3829 3996  
Fax: +61 (0)7 3829 3997  
sales@kinesis-australia.com.au  
Web: kinesis-australia.com.au

### INDIA

#### Kinesis India Pvt. Ltd

Tel: +91-22-661 41 668  
Fax: +91-22-6668 1600  
info@kinesisindia.net  
Web: kinesisindia.net

Kinesis recognise all trademarks and registered trademarks. E&OE



**Kinesis Ltd** [kinesis.co.uk](http://kinesis.co.uk)  
**Kinesis Inc** [kinesis-usa.com](http://kinesis-usa.com)  
**Kinesis GmbH (formerly Abimed)** [kinesisgmbh.de](http://kinesisgmbh.de)  
**Kinesis Australia Pty Ltd** [kinesis-australia.com.au](http://kinesis-australia.com.au)  
**Kinesis India Pvt. Ltd** [kinesisindia.net](http://kinesisindia.net)

 **Kinesis**  
Scientific Experts