

Soxhlet Extraction Thimbles

- Manufactured from pure cellulose, neutral cotton
- One-piece moulded
- Fat and impurity free
- Good mechanical strength
- Comprehensive range of sizes
- Suitable for Soxhlet extraction of most solid sample types



Soxhlet Extraction Thimbles

KX Soxhlet Extraction Thimbles are manufactured to meet the demands of today's analytical applications.

Manufactured from pure cellulose, KX Soxhlet Extraction Thimbles are used in a variety of applications where the original sample is solid in nature; common applications include the extraction of fats or pesticides in food and soil materials.

Kinesis offers an extensive range of sizes, ranging from 37 x 21 to 250 x 80mm.

The one piece moulded design offers excellent mechanical strength and rigidity, ensuring retention of the sample in the thimble throughout the extraction process.

Ordering Information

The following table provides details of the full range of KX Soxhlet Extraction Thimbles.

Part Number	Description*	Pack Size
SET-018-021-037	KX Soxhlet Thimble, 18 x 21 x 37mm	25
SET-022-025-060	KX Soxhlet Thimble, 22 x 25 x 60mm	25
SET-022-025-080	KX Soxhlet Thimble, 22 x 25 x 80mm	25
SET-025-029-080	KX Soxhlet Thimble, 25 x 29 x 80mm	25
SET-026-030-060	KX Soxhlet Thimble, 26 x 30 x 60mm	25
SET-026-030-080	KX Soxhlet Thimble, 26 x 30 x 80mm	25
SET-027-031-080	KX Soxhlet Thimble, 27 x 31 x 80mm	25
SET-030-034-060	KX Soxhlet Thimble, 30 x 34 x 60mm	25
SET-030-034-080	KX Soxhlet Thimble, 30 x 34 x 80mm	25
SET-030-034-100	KX Soxhlet Thimble, 30 x 34 x 100mm	25
SET-033-038-080	KX Soxhlet Thimble, 33 x 38 x 80mm	25
SET-033-038-130	KX Soxhlet Thimble, 33 x 38 x 130mm	25
SET-037-042-130	KX Soxhlet Thimble, 37 x 42 x 130mm	25
SET-037-042-150	KX Soxhlet Thimble, 37 x 42 x 150mm	25
SET-041-046-123	KX Soxhlet Thimble, 41 x 46 x 123mm	25
SET-041-046-150	KX Soxhlet Thimble, 41 x 46 x 150mm	25
SET-046-051-130	KX Soxhlet Thimble, 46 x 51 x 130mm	25
SET-046-051-165	KX Soxhlet Thimble, 46 x 51 x 165mm	25
SET-054-059-170	KX Soxhlet Thimble, 54 x 59 x 170mm	25
SET-058-062-170	KX Soxhlet Thimble, 58 x 62 x 170mm	25
SET-064-070-240	KX Soxhlet Thimble, 64 x 70 x 240mm	25
SET-070-077-200	KX Soxhlet Thimble, 70 x 77 x 200mm	25
SET-070-077-240	KX Soxhlet Thimble, 70 x 77 x 240mm	25
SET-080-087-200	KX Soxhlet Thimble, 80 x 87 x 200mm	25
SET-080-087-250	KX Soxhlet Thimble, 80 x 87 x 250mm	25

*The (nominal) dimensions given are the internal diameter, external diameter and length respectively.

For additional information on KX Soxhlet Extraction Thimbles or other Kinesis Sample Preparation Products visit kinesis-group.com.

KLTD-1500-0615



Kinesis recognise all trademarks and registered trademarks. E&OE


Kinesis
 Scientific Experts

sales@kinesis-group.com
kinesis-group.com